



SPECIAL
TOMATO

Multi-Positioning Seat

Special Tomato Multi-Positioning Seat– Adapt and Adjust the System as your child's needs change!



The Special Tomato Multi-Positioning Seat has been designed to offer mobile, comfortable and proper positioning support at home, school, therapy or day care.

MPS Base was designed specifically for the MPS Seat and allows for 10° - 25° of Tilt-In-Space, this accommodates almost every range of positioning needs. The MPS Mobile Floor Kit has proved valuable at school, home and therapy!

The Special Tomato MPS Mobile Floor Kit is a soft, lightweight and durable seating surface with cushions that allow for customized support.

The cushions move up and down on a hook and loop track allowing for growth of your child. The cushions are made from a latex-free material that is peel and tear resistant and impermeable to fluids.

The cushions are supported by the outer shell that maintains position. Included with the Special Tomato MPS Mobile Floor Kit is Hip-Flex adjustment that provides 30° of pivot. Bring home the seat your child needs to feel part of class and to have fun with friends!

The Special Tomato MPS has been successfully tested and Approved for aircraft travel (FAA approved) this includes all systems ordered April 1, 2006 and later.

MPS Standard Features

Seat Shell

Adjustable standard seat and back cushions

Adjustable footrest

LATCH System

Adjustable 5-Point positioning H-harness

Wooden base with handle

Four swivel castors with rear locking brakes

*regular or lateral head support must be selected-not inc as std item



Choice of two sizes Small or Large and three modern colours, **Chocolate** with **Choc** Seating, **Chocolate** with **Aqua Blue** Seating or **Chocolate** with **Lilac** Seating

MPS Optional Accessories



MPS Tray

The tray is depth adjustable and offers 60° of tilt, this will provide the best tray position for any activity like eating, drawing and reading. **The tray may not be used during transportation.** It can be stored on the back of the Mobile Floor Kit when not in use.

Available in **Chocolate** only.



Regular Head Support

The Regular Head Support is made from the same material as the MPS Seat and offers the same Anti-microbial protection and peel and tear resistance.

If your child requires more lateral support, the Lateral Head Support option would be a better fit. Available in **Chocolate**, **Aqua Blue** or **Lilac**.



Extended Seat Cushion

As your child begins to outgrow the Large MPS, adding the optional Extended Seat Cushion can extend the life of your MPS. The ideal time to purchase the Extended Seat Cushion is when your child's thigh

depth and hip width have outgrown the standard seat cushion.

You can extend the life of your MPS by 7 years with the addition of the Extended Seat Cushion.

Available in **Chocolate**, **Aqua Blue** or **Lilac**.



Lateral Head Support

The Lateral Head Support offers extra support for your child while sitting in the MPS Seating System. It is made from the same material as the MPS Seat and offers the same Anti-microbial protection and peel and tear

resistance. If your child requires less lateral support, the Standard Head Support option would be a better fit.

Available in **Chocolate**, **Aqua Blue** or **Lilac**.

Special Tomato MPS Wooden Mobile Base Sizing Chart

	Fits Small and Large MPS Sitter
Maximum User Weight	68 kg
Width for Small configuration	56 cm
Width for Large configuration	53 cm
Length	72 cm
Height	66 cm
Weight of Mobile base	13.6 kg
Seat Tilt-In-Space	10° - 25°



Special Tomato MPS Sitter Sizing Chart

	Small	Small	Large	Large
	Regular Seat	Extended Seat	Regular Seat	Extended Seat
Maximum User Weight	36.3 kg	36.3 kg	68 kg	68kg
Target Height	76-127 cm	76-127 cm	122-152 cm	122-152 cm
Target Age	2 to 5	5 to 8	8 to 11	11 to 14
Seat Depth	20.3 cm	23.75 cm	30.5 cm	39.4 cm
Seat Width	40.6 cm	40.6 cm	30.5 cm	33 cm
Seat to footrest	15.2-40.6 cm	15.2-40.6 cm	45.7 cm	45.7 cm
Top of back cushion	35.6-48.3 cm	35.6-48.3 cm	45.7 cm	45.7 cm
Top of head cushion	53.3-40.6 cm	33.0-48.3 cm	66-76.2 cm	66-76.2 cm
H-Harness height adjustment	33-48.3 cm	33-48.3 cm	53.3 cm	56.3 cm
Seat back angle (on MPS base)	-15° to + 15°	15° to + 15°	15° to + 15°	15° to + 15°
Outside Seat Width	40.6 cm	40.6 cm	47 cm	47 cm
Floor to seat height (on MPS base)	35.6 cm	35.6 cm	35.6 cm	35.6 cm
Activity Tray - min work surface	55.9w x 33d cm	55.9w x 33d cm	55.9w x 33d cm	55.9w x 33d cm



Sydney Head Office - Ph: (02) 9838 8869 Email: info@dejay.com.au

Victorian Office - Ph: (03) 9558 9669 Email: jackie@dejay.com.au - Queensland Office - Ph: (07) 3863 2589 Email: malcolm@dejay.com.au

Importers & Manufacturers of Rehabilitation Mobility Equipment